



# Padzil Dummy

Silicone Distal Pads, for Lamination Only



**ALPS Padsil** is made from the non-allergenic properties of silicone. ALPS Padsils provide more cushioning, softness and durability for your patient due to the unique compression set characteristics of silicone, these pads will not take a set even with prolonged hard use. Padsils can be used in practically every type of socket, even inside silicone suction sockets and liners. Padsils are machinable and can be shaped to fit with standard tools. **ALPS Padsil (PDS-xx) is available in 5 sizes S, M, L, XL, and XXL.** Padsil Dummy instructions follow the sizing information listed below.

### Sizing Instructions

**Step 1:**

Measure the circumference in cm at the beginning of the cylindrical section approximately 6 cm from the bottom of the distal end. This is the measured size.

**Step 2:**

Find the measured Size from Step 1 in the following chart. Read across to the find the recommended size.

<b>Limb Circumference (cm)</b>	<b>PadSil</b>
<b>14-19 cm</b>	<b>S</b>
<b>20-23 cm</b>	<b>M</b>
<b>24-27 cm</b>	<b>L</b>
<b>28-31 cm</b>	<b>XL</b>
<b>32-35 cm</b>	<b>XXL</b>


### ALPS Blue Padsil Dummies

Once the proper size has been determined, ALPS Blue Padsil Dummies are to be used when laminating or pulling plastics over the positive model to create a socket. Simply place the Blue Padsil Dummy on the distal end of the cast, and laminate or pull your plastic in the normal manner. Once the socket

is complete, pull the Blue Padsil Dummy out of the socket and replace it with a standard Padsil. These Padsils are blue so they are not mistaken for the standard Padsil, and they are made of a much stiffer durometer silicone so there is no deformation of the Padsil when creating the socket. The Blue Padsil Dummies can be used repeatedly.

### **Cleaning the Padsil Dummies**

Wash the padsil by hand with a mild soap; such as Alps Cleanser (PD 595). Allow them to dry before the next use. Be sure **not** to use any antibacterial or deodorant soap on the padsil for this may damage the product.



### **ALPS Satisfaction Guarantee**

ALPS offers a no hassle return policy within **30 days** from date of purchase. If you are unsure of the proper size of Padsil to order, you can order two sizes from ALPS Customer Service. Once received, decide which size fits best. You can then call ALPS Customer Service for a Return Authorization Number and return the improper size for full credit of the product cost within 30 days from date of purchase (excludes freight charges).

## Warranty

ALPS Padsils are warranted against manufacture's defects for **6 months** from the date of purchase. ALPS warrants only that its products will meet its specifications. **THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR USE, NOR ANY OTHER EXPRESS OR IMPLIED WARRANTIES.** The user's exclusive remedy and ALPS' sole liability is limited to the replacement of any product shown to be otherwise than as warranted. ALPS will not be liable for incidental or consequential damages of any kind.

To obtain an ALPS Return Authorization Number (RA#) call ALPS Customer service between the hours of 8 am and 5 pm EST and provide ALPS representative with the following information.

1. **Customer ID Number**
2. **Invoice Number**
3. **Date of Invoice**
4. **Nature of the return**

The RA# must be displayed on the exterior of the returned item box or it will be refused at the dock.



CE

EC	REP
----	-----

**CEpartner4U**  
Esdoornlaan 13, 3951DB Maarn, The Netherlands  
[www.cepartner4u.com](http://www.cepartner4u.com)

900028 Rev G, CAN# 6652  
December 2016  
PDSB

**ALPS South LLC**  
**2895 42<sup>nd</sup> Ave. N.**  
**St. Petersburg, FL. 33714**  
**USA**  
**Tel: 1-727-528-8566**  
**Fax: 1-727-528-8862**  
[Info@easyliner.com](mailto:Info@easyliner.com)  
[www.easyliner.com](http://www.easyliner.com)