

# SKIN RELIEVER

## ENCP

Gel Sheath - No Fabric or Distal Attachment



### DESCRIPTION

The ALPS™ Skin Reliever, made with EasyGel, eliminates shear forces, friction, and prevents abrasions on the skin of the stump. It's tapered structure accommodates volume fluctuations over the life of the prosthesis. It provides superior elongation for very little compression against the skin, reducing abrasions. ALPS Skin Reliever can be worn under silicone or gel, and can be used to replace socks.

### TOP FEATURES & BENEFITS

- ✓ Features ALPS EasyGel
- ✓ Replaces socks
- ✓ Reduces shearing and abrasions
- ✓ High elasticity
- ✓ May soften and soothe skin

#### PRODUCT NUMBER

ENCP-XX

Skin Reliever, Tapered from 3mm to 1.5mm


\*XX indicates size. See sizing chart and illustrations on reverse for further details.



#makinglivesbetter


MAKING LIVES BETTER

# CROSS SECTION DIAGRAM



## HOW TO SIZE

Measurements taken  
6cm from the distal end



## SIZING CHART

### MEASURED SIZE (CM)

10-19 cm	ENCP16
20-27 cm	ENCP20
28-37 cm	ENCP28



T: 727.528.8566 | 800.574.5426 | F: 727.528.8862

[www.easyliner.com](http://www.easyliner.com) | [info@easyliner.com](mailto:info@easyliner.com)

2895 42nd Avenue North St. Petersburg, Florida 33714

