

PADSIL™

PDS



Silicone Distal Pad

DESCRIPTION

The Padsil™ was developed to increase patient comfort in combination with ALPS® liners. The more active user will appreciate the additional cushioning, softness, and durability. Users who experience volume changes will also appreciate the Padsil's ability to lift the limb back into the proper position. They can be placed at the bottom of the socket or inside a pin locking liner.

TOP FEATURES & BENEFITS

- ✔ Soft medical grade silicone
- ✔ Helps maintain total contact
- ✔ Blue dummy padsils are available for laminating
- ✔ Keeps its shape over time
- ✔ Easily modified

PRODUCT NUMBER

PDS-XX Silicone Distal Pad

PDSB-XX Blue Dummy Distal Pad

*XX indicates size. See sizing chart and illustrations on reverse for further details.


FITTING DIAGRAM

! Can be used with any system

Can be placed inside Pin Locking Liners



Can be placed inside the socket or outside Cushion Liners



PADSIL DUROMETER


PART	SIZE	DUROMETER
PADSIL(PDS)	S-XXL	35-55
PADSIL DUMMY (PDSB)	S-XXL	65-85

SIZING CHART

MEASURED SIZE (CM)	
14-19 cm	PDS-S
20-23 cm	PDS-M
24-27 cm	PDS-L
28-31 cm	PDS-XL
32-35 cm	PDS-XXL

HOW TO SIZE

Measurements taken 6cm from the distal end



T: 727.528.8566 | 800.574.5426 | F: 727.528.8862
www.easyliner.com | info@easyliner.com
 2895 42nd Avenue North St. Petersburg, Florida 33714

